

HP PTA Balloon Dilatation Catheter



Size Matrix of HP PTA Balloon Dilatation Catheter (0.0350" OTW DLT)

Diameter (mm)	Length (mm)									
	20	30	40	60	80	100	120	150	180	200
3.0	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
4.0	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
5.0	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
6.0	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
7.0	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
8.0	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
9.0	✓	✓	✓	✓	✓	-	-	-	-	-
10.0	✓	✓	✓	✓	✓	-	-	-	-	-
12.0	✓	✓	✓	✓	✓	-	-	-	-	-

Product Features

- Wide variety of diameters and lengths up to 200mm.
- Low profile enables tight tracking to the wire.
- 6-wings balloon fluting ensures the balloon has a minimum profile to cross the stenosis easier.
- Balloon material and design allow for shape retention at rated burst pressure.
- Higher pressure when needed (up to 24 atm).
- Hydrophilic shaft coating to increase trackability.